




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine temperature sensor, PT1000-A, stainless steel

**€79.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Temperature sensor as an accessory for magnetic stirrers with heating.  
The temperature range is -50 to 800 °C.  
In stainless steel version.

## Additional Information

---

No.	D-6014
Manufacturer (Brand)	neoLabLine
EAN	4058072017934
Transport temperature	Room temperature
diameter	3 mm
Material	Stainless steel
Temp MAX	800 °C
Temp MIN	-50 °C
Length	230 mm
Measuring accuracy °C	1 °C
Type temperature control	PT1000

