



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine temperature sensor for magnetic stirrer, glass-coated

€179.40
plus VAT &
Shipping

Product Images



Description

Temperature sensor as an accessory for magnetic stirrers with heating.
The temperature range is -50 to 800 °C.
In the glass-coated version.

Additional Information

No.	D-6041
Manufacturer (Brand)	neoLabLine
EAN	4058072158200
Transport temperature	Room temperature
diameter	6 mm
Material	Borosilicate glass
Temp MAX	800 °C
Temp MIN	-50 °C
Length	230 mm
Measuring accuracy °C	1 °C
Type temperature control	PT1000

