



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN work protection hood over keypad and housing standard, retrofittable, for models with weighing plate size Ø 105 mm (fig. similar)

€44.00
plus VAT &
Shipping

Product Images



Description

- Dust cover over keypad and housing as standard, can be retrofitted, for models with weighing plate size Ø 105 mm (Fig. similar)

Additional Information

No.	KP-5339
Manufacturer (Brand)	KERN
EAN	4058072148676
Transport temperature	Room temperature
Color	transparent (milky)
Material	Polyvinyl chloride (PVC)
Manufacturer	CORE

