



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag Extension cable MAXdrive Accessories for MAXdrive/FABdrive

€131.00
plus VAT &
Shipping

Product Images



Description

Extension cable for MAXdrive and FABdrive, length 3m, temperature resistant from -10°C to +50°C. Shipping weight (gross): approx. 0,27 kg

Additional Information

No.	ZU-0045
Manufacturer (Brand)	2mag
EAN	4058072444921
Gross weight	0.27
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.21000 kg
Country of origin	Germany
Packing width	0.15000 m
Packing height	0.10000 m
Packing depth	0.15000 m
Packaging volume	0.00225000 cc
Customs tariff number	85442000

