



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag extension cable MIXdrive Accessories for MIXdrive

€97.00
plus VAT &
Shipping

Product Images



Description

Extension cable for MIXdrive, length 3m, temperature resistant from -10°C to +50°C. Shipping weight (gross): approx. 0.15 kg

Additional Information

| | |
|-----------------------|------------------|
| No. | ZU-0043 |
| Manufacturer (Brand) | 2mag |
| EAN | 4058072444891 |
| Gross weight | 0.15 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| Net weight | 0.09000 kg |
| Country of origin | Germany |
| Packing width | 0.15000 m |
| Packing height | 0.10000 m |
| Packing depth | 0.15000 m |
| Packaging volume | 0.00225000 cc |
| Customs tariff number | 85442000 |

