




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag bioMIXcontrol 4MS Control unit for bioMIXdrive 1

€2,675.00
plus VAT &
Shipping

Product Images



Description

Powerful control unit, for 4 stirrer drives bioMIXdrive 1, developed for the demanding daily laboratory routine and uncompromising heat-free continuous operation. Control unit for the connection of up to 4 stirrer drives bioMIXdrive 1, for each stirrer drive individually and independently adjustable speeds (master-slave technology), as well as ON/OFF function, synchronous operation of all drives at the push of a button, extended speed range from 5 to 250 rpm, heat-free 6 W per stirring point, automatic power setting for heat-free long-term operation, jerk-free and uniform stirring even at the lowest speeds and reduced power for heat-free continuous operation e.g. in CO2 incubators. Digital display for speed setting, SoftStart for reliable and safe acceleration of the stirring staff, QuickSet for fast setting of the initial or maximum speed, saving of the last functional state when switched off, durable stainless steel housing, easy to clean. 3 years warranty on all material and manufacturing defects, developed and made in Germany. Speed range: 5 - 250 rpm (4x individually adjustable) Connection sockets: 4 Stirring power (max.): 6 W / stirring point, automatic Power setting: none Material housing: stainless steel Perm. operating conditions: 0 up to +40 °C (at 80% humidity) Perm. storage conditions: -40 °C to +40 °C (at 80% humidity) Perm. storage conditions: -40 °C to +40 °C (at 80% humidity) Storage conditions: -40 °C to +70 °C, 10 - 95 %, 500 - 1060 hPa Connected load: 100-240 V/50-60 Hz/1.5 A Output voltage: 48 VDC Protection class: IP20 Dimensions (WxDxH): 200 x 155 x 38 mm Shipping weight (gross): approx. 2.0 kg

Additional Information

No.	ZU-0087
Manufacturer (Brand)	2mag
VGKL number	101859864
EAN	4058072445232
Gross weight	2
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.70000 kg
Country of origin	Germany
Packing width	0.29000 m
Packing height	0.15000 m
Packing depth	0.24000 m
Packaging volume	0.01044000 ccm
Customs tariff number	90328900

