



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag distriBOX 8 Accessories for MIXdrive distributor box

€348.00
plus VAT &
Shipping

Product Images



Description

8-fold distributor for synchronous operation of 8 stirrer drives MIXdrive on one control unit MIXcontrol 20 or MIXcontrol 40, identical operating conditions, synchronous speed, same stirring performance, space-saving, potted stainless steel housing. Dimensions: 48 x 92 x 25 mm Shipping weight (gross): approx. 0.27 kg

Additional Information

No.	ZU-0051
Manufacturer (Brand)	2mag
VGKL number	101854908
EAN	4058072444976
Gross weight	0.27
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.21000 kg
Country of origin	Germany
Packing width	0.15000 m
Packing height	0.10000 m
Packing depth	0.15000 m
Packaging volume	0.00225000 cc
Customs tariff number	90328900

