



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## 2mag mixBATH accessories for MIXdrive MIXdrive1/1 eco/1 XS//6/15/60 stainless steel

**€794.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Stainless steel bath, suitable for MIXdrive 6/15/60 as well as MIXdrive 1/1 eco/1 XS, durable, stable and easy to clean, suitable for immersion circulators (heating/cooling), easy to clean due to gap-free welded edges, integrated drain cock for easy and safe draining, electropolished, all-round edge designed as carrying handle, exact positioning of the stirrer for underbath operation, stirring point marking in the bath, left/right stop for optimum thermostat positioning (circulation pump) when operating with MIXdrive 6 or MIXdrive 15. Inner dimensions (WxDxH): 277 x 537 x 175 mm Outer dimensions (WxDxH): 318 x 600 x 190 mm Shipping weight (gross): approx. 8.0 kg

## Additional Information

---

No.	ZU-0031
Manufacturer (Brand)	2mag
VGKL number	101854950
EAN	4058072444778
Gross weight	8
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	6.00000 kg
Country of origin	Germany
Packing width	0.64000 m
Packing height	0.27000 m
Packing depth	0.40000 m
Packaging volume	0.06912000 cc
Customs tariff number	73102990

