



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag mixBATHcover for MIXdrive 15 Accessories for MIXdrive 15, cover mixBATH

€321.00
plus VAT &
Shipping

Product Images

Description

Chemically resistant, transparent polycarbonate lid for mixBATH or mixBATH PC, to minimize evaporation effects. Different cut-outs available for positioning the thermostat in bath operation of MIXdrive 15 and MIXdrive 6. Material: polycarbonate Dimensions (WxDxH): 315 x 575 x 8 mm Shipping weight (gross): approx. 2.6 kg

Additional Information

No.	ZU-0029
Manufacturer (Brand)	2mag
VGKL number	101854815
EAN	4058072444754
Gross weight	2.6
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.70000 kg
Country of origin	Germany
Packing width	0.63000 m
Packing height	0.33000 m
Packing depth	0.32000 m
Packaging volume	0.06652800 cc
Customs tariff number	39239000

