



Bandelin D 1050 CS lid

Product Images



Description

A suitable stainless steel lid for an ultrasonic bath protects the bath liquid from external contamination and prevents odor nuisance caused by the added cleaning or disinfecting preparations. In addition, the noise emission during operation of the ultrasonic cleaner, which may be perceived as annoying, is reduced. The lids are designed in such a way that any condensation water produced is directed back into the tank.

No.	BL-0505
Manufacturer (Brand)	Bandelin
EAN	4058072504045
Country of origin	Deutschland
Region of origin	Berlin
Customs tariff number	84799070



Additional Information