



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin ED 9 Injection nozzle holder

€86.00
plus VAT &
Shipping

Product Images



Description

ED 9 injection nozzle holder made of stainless steel. For mounting in the SONOREX ultrasonic cleaner. Accommodates max. 32 injection nozzles: 8 holes diameter: 9.5 mm each 8 holes diameter: 5.5 mm each 8 holes diameter: 4.5 mm each 8 holes diameter: 3.5 mm each

Additional Information

No.	BL-0021
Dual use indicator	false
Manufacturer (Brand)	Bandelin
EAN	4058072503673
Transport temperature	Room temperature
Country of origin	Germany
Region of origin	Berlin
Packing width	0.23000 m
Packing height	0.06000 m
Wide	70 mm
Packing depth	0.07000 m
Packaging volume	0.00098600 cc
Customs tariff number	84799070
Height	60 mm
Length	231 mm

