




## Bandelin HG 40, holding frame

### Product Images

---



### Description

---

A holding frame is used for the correct and variable positioning of the ultrasonic transducers. - Suitable for direct and indirect sonication - Application possible with and without noise protection box LS 40 - Commercially available sleeves can be mounted on stand rod (d=16 mm) - Material base and rod: stainless steel

# Additional Information

|                       |                |
|-----------------------|----------------|
| No.                   | BL-0495        |
| Dual use indicator    | false          |
| Manufacturer (Brand)  | Bandelin       |
| VGKL number           | 100733681      |
| EAN                   | 4058072506100  |
| Country of origin     | Deutschland    |
| Region of origin      | Berlin         |
| Packaging volume      | 0.01045000 ccm |
| Customs tariff number | 84799070       |

