



Bandelin KD 0, SIEBKORB

Product Images



Description

Stainless steel, sieve mesh with very fine mesh for insertion into insertion vessels or for cleaning very fine parts directly in the ultrasonic bath.

Additional Information

No.	BL-0255
Dual use indicator	false
Manufacturer (Brand)	Bandelin
VGKL number	960736202
EAN	4058072506155
Country of origin	Deutschland
Region of origin	Berlin
Packing width	0.08000 m
Packing height	0.11000 m
Packing depth	0.08000 m
Wide	82 mm
Customs tariff number	84799070
Height	40 mm
Length	82 mm

