



Bandelin MS 72, MICROSPITZE

Product Images



Description

Sonotrodes transfer the mechanical vibrations into the sample, are thermostable, autoclavable and resistant to corrosive media. Sample volume, vessel diameter and desired amplitude determine device selection and sonotrode type. A high amplitude means particularly intensive sonication. The sonotrode is connected to step or booster horns. All BANDELIN sonotrodes are made of the titanium alloy TiAl6V4 (3.7165). Micro tips are sonotrodes with a particularly thin tip, which enables sonication of small volumes with a high amplitude and thus high ultrasonic effect.

Additional Information

No.	BL-0307
Reg. trademark	SONOPULS ®
Dual use indicator	false
Manufacturer (Brand)	Bandelin
VGKL number	530703062
EAN	4058072506803
diameter	2 mm
Country of origin	Deutschland
Region of origin	Berlin
Packing width	0.23000 m
Packing height	0.03000 m
Packing depth	0.07000 m
Packaging volume	0.00047300 ccm
Customs tariff number	84799070

