



Bandelin PD 04, SIEBKORB

Product Images



Description

Strainer basket with very fine mesh for insertion into insertion vessels or for cleaning very fine parts directly in the ultrasonic bath.

Additional Information

No.	BL-0008
Dual use indicator	false
Manufacturer (Brand)	Bandelin
VGKL number	960736204
EAN	4058072503567
Country of origin	Deutschland
Region of origin	Berlin
Packing width	0.07000 m
Packing height	0.10000 m
Packing depth	0.07000 m
Packaging volume	0.00046600 ccm
Wide	75 mm
Customs tariff number	39269097
Height	50 mm
Length	75 mm

