



Bandelin TS 113, Titanium sonotrode

Product Images



Description

Diameter 13 mm fits SH 100 G and SH 200 G (HD 4100 and HD 4200). Sonotrodes transfer the mechanical vibrations into the sample, are thermostable, autoclavable and resistant to corrosive media. Sample volume, vessel diameter and desired amplitude determine instrument selection and sonotrode type. A high amplitude means particularly intensive sonication. The sonotrode is connected to step or booster horns. All BANDELIN sonotrodes are made of the titanium alloy TiAl6V4 (3.7165).

No.	BL-0274
Dual use indicator	false
Manufacturer (Brand)	Bandelin
EAN	4058072506346
diameter	13 mm
Country of origin	Deutschland
Region of origin	Berlin
Packing width	0.25000 m
Packing height	0.03000 m
Packing depth	0.03000 m
Customs tariff number	84799070



Additional Information