




Bandelin TT 13, TITANTELLER

Product Images



Description

Diameter 13 mm fits SH 70 G (HD 2070.2) and SH 213 G (HD 2200.2). A titanium disc represents a greatly shortened sonotrode. The ultrasonic emitting surface is subject to wear by cavitation erosion in any design. However, the replacement of a titanium disk can be done at a much lower cost than for long sonotrodes due to the low material cost. The titanium disks are used for medium volume sonication.

No.	BL-0310
Reg. trademark	SONOPULS ®
Dual use indicator	false
Manufacturer (Brand)	Bandelin
VGKL number	530705215
EAN	4058072506834
diameter	13 mm
Country of origin	Deutschland
Region of origin	Berlin
Packing width	0.18000 m
Packing height	0.03000 m
Packing depth	0.10000 m
Customs tariff number	84799070



Additional Information