



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



SmartRack® storage box insert - for centrifuge tubes

€59.99
plus VAT &
Shipping

Product Images



Description

The module can be inserted precisely into the storage box and enables easy storage and fixation of centrifuge tubes.

Additional Information

No.	BS-0010
Manufacturer (Brand)	SmartRack
EAN	4066292683754
Transport temperature	Room temperature
suitable for	Centrifuge tube

