



## Better Basics SmartCase

### Product Images



### Description

"

The **SmartCase** is designed to transport and store the SmartRack®, a modular organization system for laboratory workstations. Its aluminum case safely protects a SmartRack® inside and its contents from light and other environmental influences.

The SmartCase consists of a lower part with four castors into which a SmartRack® can be placed. The upper part of the SmartCase has a cushion on the inside which serves to fix the SmartRack® inside the SmartCase as soon as the SmartCase is closed by the tension locks. This prevents the SmartRack® from slipping or tipping over.

The SmartCase has a self-acting valve system via which the atmosphere inside the SmartCase can be changed by an external gas supply. This allows, for example, the temporary creation of a low-oxygen atmosphere in the SmartCase, which creates an additional safety level when transporting oxygen-sensitive, highly flammable or spontaneously combustible substances in air.

**The SmartCase comprises the following components:**

A bottom section with four casters, four tension locks and an automatic gas release valve system

A top section with a handle, four touchdown latches and a self-actuating gas inlet valve system

A connector for connecting gas hoses of various diameters to the SmartCase.

"

# Additional Information

No.	BS-0050
Manufacturer (Brand)	Better Basics
EAN	4066292684157

