



Bühler spring clamp 500 ml

Product Images



Description

For mounting on the universal trays. Sizes matched to Erlenmeyer flasks, but can also be used for round-bottomed flasks, beakers, etc.

Additional Information

No.	NB-0177
Dual use indicator	false
Manufacturer (Brand)	Bühler
EAN	4058072623890
Gross weight	0.5
Net weight	0.11000 kg
Material	Edelstahl
Country of origin	Deutschland
Region of origin	Baden-Württemberg
Customs tariff number	84799070

