



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Wide Neck Bottle, 1.000 mL, for Rotor S-4xUniversal, 2 pcs

€222.00
plus VAT &
Shipping

Product Images



Description

Wide neck bottle, 1,000 mL, for rotor S-4xUniversal, set of 2 - Capacity: 1,000 mL - Shape: round, wide neck, flat bottom - Size (diameter x height): 118 x 133 mm - Life: at least 2 years or 250 cycles à 30 min at 4 °C - 25 °C - Autoclavable: at least 50 cycles at 125 °C - Material: polypropylene (PP) - Max. RZB: 4,415 x g (Centrifuge 5910 R with rotor S-4xUniversal) - Set of 2 pieces - Adapter required (Catalog No. 5910 756 003)

Additional Information

No.	VB-2547
Manufacturer (Brand)	Eppendorf
EAN	4058072425562
Gross weight	0.55
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.49500 kg
Country of origin	Germany
Customs tariff number	39269097

