




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Adapter, for bottle 500 mL Corning®, 2 pieces

**€313.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Adapter, for bottle 500 mL Corning®, for rotor S-4-104, rotor S-4x1000 round cup and rotor S-4x750, 2 pieces

## Additional Information

---

No.	VB-2264
Manufacturer (Brand)	Eppendorf
EAN	4058072422974
Gross weight	0.5
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.45000 kg
Country of origin	Germany
Packing width	0.19000 m
Packing height	0.14000 m
Packing depth	0.19000 m
Packaging volume	6000.000000 ccm
Customs tariff number	84219100

