



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf rotor lid FA-45-16-17, aerosol-tight, for Centrifuge 5430/5430R

€252.00
plus VAT &
Shipping

Product Images



Description

Rotor lid for rotor FA-45-16-17, aerosol-tight, aluminum, with QuickLock, for Centrifuge 5430/5430R

Additional Information

No.	VB-1940
Manufacturer (Brand)	Eppendorf
EAN	4058072419851
Gross weight	0.35
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.31500 kg
Country of origin	Germany
Packing width	0.18000 m
Packing height	0.06000 m
Packing depth	0.27000 m
Packaging volume	3000.000000 ccm
Customs tariff number	84219100

