




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf replacement clamp, for adapter for 250 mL rectangular beakers, 2 pcs

€9.10
plus VAT &
Shipping

Product Images



Description

Replacement clamp, for adapter for 250 mL rectangular beakers, 2 pcs

Additional Information

No.	VB-2183
Manufacturer (Brand)	Eppendorf
EAN	4058072422202
Gross weight	0.044
HSNumber	84219100
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.04000 kg
Country of origin	Germany
Packing width	0.25000 m
Packing height	0.02000 m
Packing depth	0.15000 m
Packaging volume	1000.000000 ccm
Customs tariff number	84219100

