



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf thread adapter for bottles, PTFE, from GL 45 to S 40

**€29.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Thread adapter for bottles, PTFE, from GL 45 to S 40

## Additional Information

---

No.	VB-1584
Manufacturer (Brand)	Eppendorf
EAN	4058072416874
Gross weight	0.023
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.02100 kg
Country of origin	Germany
Packing width	0.14000 m
Packing height	0.03000 m
Packing depth	0.08000 m
Packaging volume	1000.000000 ccm
Customs tariff number	84799070

