



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Tube Holder 1.5/2.0 mL, for 24 × 1.5 mL or 2.0 mL reaction tubes

€110.00
plus VAT &
Shipping

Product Images



Description

Vessel holder 1.5/2.0 mL, for 24 × 1.5 mL or 2.0 mL reaction vessels - High-speed mixer (3.000 rpm) for fast and effective mixing of small volumes in all common laboratory tubes up to 50 mL - 2DMix-Control technology (planar mixing motion) for controlled mixing and reproducible mixing results in seconds - Anti-Spill technology for controlled mixing motions without lid wetting and cross-contamination - 5 different tube holders (0.5 mL, 1.5/2.0 mL, 5/15 mL, 25/50 mL, PCR-96 well plates) for high vessel flexibility - Integrated ergonomic vortex function for individual sample mixing of different vessel formats - Ergonomic operation due to minimal hand-arm vibration transmission during vortexing (certificate available) (PhysioCare Concept®) - 5 function keys with predefined standard mixing parameters for easy operation - High stability - remains stable even at highest mixing speed - Extremely quiet operation to improve working environment

Additional Information

No.	VB-1791
Manufacturer (Brand)	Eppendorf
VGKL number	106934130
EAN	4058072418526
Gross weight	0.2
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.18000 kg
Country of origin	Germany
Packing width	0.21000 m
Packing height	0.07000 m
Packing depth	0.14000 m
Packaging volume	3000.000000 ccm
Customs tariff number	84799070

