



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Eppendorf thread adapter, for canister 5 L, ETFE, from 45 mm to 17/8 and No. 039

**€60.80**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Thread adapter, for canister 5 L, ETFE, from 45 mm to 17/8and#039;thread

## Additional Information

---

No.	VB-1582
Manufacturer (Brand)	Eppendorf
EAN	4058072416850
Gross weight	0.094
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.08500 kg
Country of origin	Germany
Packing width	0.24000 m
Packing height	0.07000 m
Packing depth	0.15000 m
Packaging volume	3000.000000 ccm
Customs tariff number	84799070

