



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf adapter plate for 96 × 0.2 mL PCR tubes and PCR plates 96, as insert in ther

**Price on  
request**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Adapter plate for 96 × 0.2 mL PCR tubes and PCR plates 96, as insert in thermoblock for MTPs – A wide range of thermoblocks for reaction tubes from 0,2 mL to 50 mL as well as plates and glass slides – Thermal insulation of the thermoblock for ergonomic working (no direct contact with hot surfaces – no burning of fingers) – The principle of dry incubation eliminates the risk of contamination that often occurs in classical water bath incubations – Easy exchange of thermoblocks – The reaction vessels are easily accessible – Extremely quiet operation and optimal temperature transfer due to precise fit and complete enclosure of the vessels in the thermoblocks

## Additional Information

---

No.	VB-1801
Manufacturer (Brand)	Eppendorf
VGKL number	466953639
EAN	4058072175245
Gross weight	0.5
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.45000 kg
Country of origin	Germany
Packing width	0.10000 m
Packing height	0.15000 m
Packing depth	0.06000 m
Packaging volume	1000.000000 ccm
Customs tariff number	84799070

