



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Tube Rack, 12 places, for 5.0 mL and 15 mL tubes, 2pcs

**€50.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Eppendorf Tube Rack, 12 places, for 5.0 mL and 15 mL tubes, polypropylene, numbered places, autoclavable, 2 pieces - Racks are stackable even with different formats - for space-saving storage in the laboratory when not in use - Non-slip silicone feet and a footprint that prevents tipping and ensures safe handling - Wave-like design: The rows of vessels are arranged at different heights, so that the reaction vessels are easily accessible and clearly visible - Made of polypropylene (PP) for high chemical and physical resistance. - Particularly space-saving due to optimal arrangement of the vessels in relation to the footprint with unrestricted handling of the vessels. - Particularly flexible: space for 5.0 mL and 15 mL tubes - Numbered tube positions for quick and easy orientation and error-free handling

## Additional Information

---

No.	VB-0302
Manufacturer (Brand)	Eppendorf
EAN	4058072405632
Gross weight	0.335
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.00000 kg
Material	Polypropylene (PP)
Country of origin	Germany
Packing width	0.25000 m
Packing height	0.12000 m
Packing depth	0.14000 m
Packaging volume	5000.000000 ccm
Customs tariff number	39269097

