



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Tube Rack, 12 places, 6 for 5.0 mL and 15 mL tubes and 6 for 25 mL u

**€50.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Eppendorf Tube Rack, 12 places, 6 for 5.0 mL and 15 mL tubes and 6 for 25 mL and 50 mL tubes, polypropylene, numbered places, autoclavable, 2 pieces – Racks are stackable even with different formats – For space-saving storage in the lab when not in use – Non-slip silicone feet and a footprint that prevents tipping and ensures safe handling – Wave-like design: The rows of vessels are arranged at different heights, so that the reaction vessels are easily accessible and clearly visible – Made of polypropylene (PP) for high chemical and physical resistance. – Particularly space-saving due to optimal arrangement of the vessels in relation to the footprint with unrestricted handling of the vessels. – Particularly flexible: vessels of different sizes fit into the 1.5/2.0 mL, 5.0/15.0 mL, 5.0/15/50 mL, and the cryovessel racks – Numbered vessel locations for quick and easy orientation and error-free handling

## Additional Information

---

No.	VB-0303
Manufacturer (Brand)	Eppendorf
VGKL number	106919835
EAN	4058072187316
Gross weight	0.363
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.00000 kg
Material	Polypropylene (PP)
Country of origin	Germany
Packing width	0.25000 m
Packing height	0.12000 m
Packing depth	0.14000 m
Packaging volume	5000.000000 ccm
Customs tariff number	39269097

