



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Tube Rack, 36 places, for cryo tubes, 2pcs

€50.20
plus VAT &
Shipping

Product Images



Description

Eppendorf Tube Rack, 36 places, for cryo tubes, polypropylene, numbered places, autoclavable, 2 pieces - Racks are stackable even with different formats - for space-saving storage in the laboratory when not in use - Non-slip silicone feet and a footprint that prevents tipping and ensures safe handling - Wave-like design: The rows of vessels are arranged at different heights, so that the reaction vessels are easily accessible and clearly visible - Made of polypropylene (PP) for high chemical and physical resistance. - Particularly space-saving due to optimal arrangement of the vessels in relation to the footprint with unrestricted handling of the vessels. - Particularly flexible: space for vessels of different sizes - Numbered vessel locations for quick and easy orientation and error-free handling

Additional Information

No.	VB-0304
Manufacturer (Brand)	Eppendorf
EAN	4058072405649
Gross weight	0.361
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.34000 kg
Country of origin	Germany
Packing width	0.25000 m
Packing height	0.12000 m
Packing depth	0.15000 m
Packaging volume	5000.000000 ccm
Customs tariff number	39269097

