



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf VisioNize® pipette manager

€1,065.00
plus VAT &
Shipping

Product Images



Description

Visionize® pipette manager, an external touch server enables communication with networked electronic pipettes. - Quickly transfer settings and other parameters to networked electronic pipettes and start pipetting right away - Let the Liquid Manager integrated into the software guide you when working with liquids that form foam or are of different viscosity or volatility - Choose from predefined Eppendorf liquid types, e.g. acetone. E.g. acetone, or create your own - Manage the pipettes in your inventory more effectively in your workgroup and improve collaboration (e.g. by sharing self-created liquid types) - All electronic Eppendorf Xplorer and Xplorer plus 4-, 6-, 8- and 12-channel pipettes are compatible. Move It® pipettes are currently not supported. - You can easily upgrade the Eppendorf Xplorer and Eppendorf Xplorer plus electronic pipettes with the Eppendorf Xplorer connect Wi-Fi module. The Wi-Fi modules have a shorter battery so that they fit optimally. - The Visionize pipette manager is based on the Visionize touch interface concept, which has been extensively tested by customers to reduce learning and adaptation times. - This is just the beginning - more features will follow.

Additional Information

No.	VB-2731
Manufacturer (Brand)	Eppendorf
EAN	4058072672195
Gross weight	1.85
HSNumber	85176200
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.66500 kg
Country of origin	Germany
Customs tariff number	85176200

