



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf thermoblocks, for 24 × 0.5 mL reaction tubes

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Thermoblocks, for 24 reaction vessels, complete with rack and cooling battery 0 °C (from the IsoTherm® system), for 24 × 0.5 mL reaction vessels - A wide range of thermoblocks for reaction vessels from 0,2 mL to 50 mL as well as plates and glass slides - Thermal insulation of the thermoblock for ergonomic working (no direct contact with hot surfaces - no burning of fingers) - The principle of dry incubation eliminates the risk of contamination that often occurs with classical water bath incubations - Easy exchange of thermoblocks - The reaction vessels are easily accessible - Extremely quiet operation and optimal temperature transfer due to precise fit and complete enclosure of the vessels in the thermoblocks

Additional Information

No.	VB-1794
Manufacturer (Brand)	Eppendorf
VGKL number	466953615
EAN	4058072169435
Gross weight	0.84
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.75600 kg
Volume MAX	0.5 ml
Country of origin	Germany
Customs tariff number	84799070

