



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf wide-neck bottle, 750 mL, for rotors S-4-104, S-4x750, S-4x1000, 2 pieces

€208.00
plus VAT &
Shipping

Product Images



Description

Wide neck bottle, 750 mL, for rotors S-4-104, S-4x750, S-4x1000, 2 pieces - Polypropylene - Centrifugation resistance up to 4,570 × g - Autoclavable

Additional Information

No.	VB-2197
Manufacturer (Brand)	Eppendorf
EAN	4058072422349
Gross weight	0.44
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.39600 kg
Country of origin	Germany
Packing width	0.20000 m
Packing height	0.15000 m
Packing depth	0.20000 m
Packaging volume	6000.000000 ccm
Customs tariff number	39269097

