



Heidolph AUTOaccurate probe

Product Images



Description

For Hei-VAP Expert Control and Hei-VAP Ultimate Control Only compatible in combination with glass set G3 or G6

Additional Information

No.	BH-1867
Manufacturer (Brand)	Heidolph
EAN	4058072779979
Gross weight	0.2
Net weight	0.20000 kg
Country of origin	Deutschland
Region of origin	Thüringen
Packing width	0.10000 m
Packing height	0.03000 m
Packing depth	0.44000 m
Packaging volume	0.00100000 ccm
Customs tariff number	90251900

