




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph Findenser™ NS24 Core, NS24 Sleeve

**€404.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Findenser™ high performance chiller in this size (400 mm) is suitable for refluxed syntheses with solvent volumes up to 1 liter. Can be mounted with standard lab clamps or with the Findenser™ mounting kit (sold separately). With the Findenser™ you can save water, space and costs when refluxing from now on!

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | BH-2191          |
| Manufacturer (Brand)  | Heidolph         |
| EAN                   | 4058072776084    |
| Gross weight          | 1.36             |
| Transport temperature | Room temperature |
| Net weight            | 1.20000 kg       |
| Country of origin     | Great Britain    |
| Customs tariff number | 84195080         |

