




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph Findenser™ NS29 core, NS24 sleeve

**€404.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Findenser™ high performance chiller in this size (400 mm) is suitable for refluxed syntheses with solvent volumes up to 1 liter. Can be mounted with standard lab clamps or with the Findenser™ mounting kit (sold separately). With the Findenser™ you can save water, space and costs when refluxing from now on!

## Additional Information

---

No.	BH-2192
Manufacturer (Brand)	Heidolph
EAN	4058072776091
Gross weight	1.38
Transport temperature	Room temperature
Net weight	1.20000 kg
Country of origin	Great Britain
Packing width	0.21000 m
Packing height	0.31000 m
Packing depth	0.21000 m
Packaging volume	0.01400000 ccm
Customs tariff number	84195080

