



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e



Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph glass set G3 XL NS 29/32

**€1,552.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Vertical coolers for all standard distillations; for maximum evaporation performance The XL coolers with a cooling surface of 2,200 cm<sup>2</sup> - achieve up to 40% higher evaporation rates. The scope of delivery includes one 1,000 ml evaporation flask and one 1,000 ml receiver flask. Standard ground joint is NS 29/32, also available with ground joint NS 24/29 on request. Including aeration cap without ground joints for easy handling and grease-free evaporation. Optional aeration valve with function for precise re-dosing available. Suitable for the automatic distillation program Dynamic AUTO accurate with Hei-VAP Expert / Ultimate Control models. Glassware sets for Soxhlet extraction and descending cooling systems are also available on request.

## Additional Information

---

No.	BH-1450
Manufacturer (Brand)	Heidolph
EAN	4058072776442
Gross weight	4.5
Transport temperature	Room temperature
Net weight	2.40000 kg
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.50000 m
Packing height	0.28000 m
Packing depth	0.51000 m
Packaging volume	0.07100000 cc
Customs tariff number	84195080

