



# Heidolph retaining bolt (SMG)

## Product Images

## Description

Retaining bolt (SMG)

## Additional Information

|                       |               |
|-----------------------|---------------|
| No.                   | BH-2494       |
| Manufacturer (Brand)  | Heidolph      |
| EAN                   | 4058072764388 |
| Gross weight          | 0.035         |
| Net weight            | 0.03500 kg    |
| Country of origin     | Deutschland   |
| Region of origin      | Bayern        |
| Customs tariff number | 73181595      |

