




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Heat-On Mutli-Well 50 ml insert

€116.00
plus VAT &
Shipping

Product Images



Description

This Heat-On insert is suitable for heating 50 ml round bottom flasks. For optimal use, it is recommended to purchase the Multi-Well holder including safety handles. The maximum temperature for continuous use is 260°C. To measure this accurately with a temperature probe, a small hole is provided on one side of the Multi-Well holder.

Additional Information

No.	BH-1370
Manufacturer (Brand)	Heidolph
VGKL number	104158630
EAN	4058072775766
Gross weight	0.265
Transport temperature	Room temperature
Net weight	0.23200 kg
Country of origin	Great Britain
Packing width	0.10000 m
Packing height	0.06000 m
Packing depth	0.11000 m
Packaging volume	0.00100000 ccm
Customs tariff number	84199085

