



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph multichannel cassette small

€261.00
plus VAT &
Shipping

Product Images



Description

Flow rates from 0.005 - 37 ml/min Suitable for tubing with 0.9 mm tubing wall thickness Available tubing diameters: 0.2 / 0.5 / 0.9 / 1.4 and 2.8 mm Special tubing piece with 2 stoppers (length: 40 cm) required for insertion into the cassette Tubing stoppers take over the tubing fixation With tubing connectors and extension tubing, it is possible to extend the tubing length by the meter Equipment: Multi-channel pump head C 4: Max. 4 x cassette small Multi-channel pump head C 12: Max. 12 x cassette small

Additional Information

No.	BH-1637
Manufacturer (Brand)	Heidolph
VGKL number	584190022
EAN	4058072778057
Gross weight	0.052
Transport temperature	Room temperature
Net weight	0.05200 kg
Country of origin	USA
Packing width	0.10000 m
Packing height	0.01000 m
Packing depth	0.14000 m
Customs tariff number	84139100

