



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Monoblock for 40 x 16 mm test tubes

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Attachment for max. 40 test tubes with 16 mm diameter. The attachment is resistant to chemicals. It is made of aluminum and is suitable for heating temperatures up to 260 °C. Together with the StarFish base plate, the MonoBlocks are part of the basic elements of the entire StarFish multi-reaction system. The MonoBlocks are attachments for several reaction vessels of the same size.

Additional Information

No.	BH-1344
Manufacturer (Brand)	Heidolph
VGKL number	104158040
EAN	4058072775506
Gross weight	0.5
Transport temperature	Room temperature
Net weight	0.10000 kg
Country of origin	Great Britain
Packing width	0.20000 m
Packing height	0.11000 m
Packing depth	0.21000 m
Packaging volume	0.00500000 ccm
Customs tariff number	84199085

