



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Monoblock for 5 x 250 ml flasks

€654.00
plus VAT &
Shipping

Product Images



Description

Attachment for max. 5 round bottom flasks 250 ml. The attachment is resistant to chemicals. It is made of aluminum and is suitable for heating temperatures up to 260 °C. Together with the StarFish base plate, the MonoBlocks are part of the basic elements of the entire StarFish multi-reaction system. The MonoBlocks are attachments for several reaction vessels of the same size.

Additional Information

No.	BH-1341
Manufacturer (Brand)	Heidolph
VGKL number	104158112
EAN	4058072775476
Gross weight	2.382
Transport temperature	Room temperature
Net weight	2.24000 kg
Country of origin	Great Britain
Packing width	0.26000 m
Packing height	0.06000 m
Packing depth	0.28000 m
Packaging volume	0.00400000 ccm
Customs tariff number	84199085

